



Microgen™ Listeria -ID

Instructions for Use

REF MID-67

MICROGEN Listeria-ID

Quick Reference

STEP 1

**SELECT A SINGLE, WELL - ISOLATED
COLONY**

STEP 2

EMULSIFY IN LISTERIA SUSPENDING BROTH

STEP 3

TRANSFER 4 DROPS TO EACH MICROWELL

STEP 4

**ADD 1 DROP OF HAEMOLYSIN REAGENT
(Well 12)**

STEP 5

INCUBATE 35 –37°C FOR 18 – 24 HOURS

STEP 6

READ AND RECORD RESULTS

STEP 7

**INTERPRET USING MICROGEN
IDENTIFICATION SYSTEM SOFTWARE**

The Microgen Listeria-ID system is intended for use by qualified laboratory personnel using aseptic technique and appropriate microbiological precautions.

The Microgen Listeria-ID system employs 12 standardised micro well substrates combined with the Microgen Identification System Software to identify members of the genus *Listeria*:

Listeria monocytogenes
Listeria welshimeri
Listeria ivanovii

Listeria innocua
Listeria grayi
Listeria seeligeri

The above organisms can be identified from selective or non-selective agar using Microgen Listeria-ID. Identification is achieved using all of the tests recommended in international standard methods for the identification of *Listeria spp.* without the need for additional confirmatory tests (1,2,3)

PRINCIPLE

Each Microgen Listeria-ID microwell strip contains 11 dehydrated substrates for the performance of carbohydrate utilisation tests and one empty well for the performance of a haemolysin reaction (4). The selection of the substrates included in the test panel is based on a combination of those substrates recommended in international standard methods (1,2,3) plus additional tests which either confirm the isolate being tested as belonging to the genus *Listeria* (Aesculin Hydrolysis, Trehalose and Arabinol Fermentation(5,6)) and/ or further enhance the differentiation of the various species comprising the genus.

Identification of isolates is achieved by recording the results visualised by a colour change after 18-24 hours incubation (there are no reagents to be added on Day 2). These results are then analysed using the Microgen Identification System Software (MID-60)

Each Microgen Listeria-ID microwell strip consists of twelve wells containing the substrates for the following 11 biochemical reactions:

		Reaction	Positive	Negative
1	Aesculin	Aesculin hydrolysis	Black	Straw colour
2	Mannitol	Fermentation of specific sugars producing acid end products changes the Bromocresol Purple indicator from purple to yellow	Yellow	Purple
3	Xylose			
4	Arabinol			
5	Ribose			
6	Rhamnose			
7	Trehalose			
8	Tagatose			
9	Glucose-1-Phosphate			
10	Methyl-D-Glucose			
11	Methyl-D-Mannose			
12	Haemolysin	Haemolysis of sheep red blood cells	Straw - Brown coloured homogeneous liquid, no carpet of red cells on the well floor	Carpet of red cells on well floor. Cells may appear red – brown in colour

Well number 12 is empty and is used for an in-well haemolysis reaction when haemolysin reagent is added to a bacterial suspension.

REAGENTS

Kit Contents (20 tests)

Holding frame for test microwell strips
Result forms
Instructions for use
20 microwell strips in individual foil pouches
20 bottles of Listeria Suspending Medium
1 bottle of Haemolysin Reagent

Additional Materials Required (not supplied in the kit)

Microgen Identification System Software (MID-60)
Sterile bacteriological loops
Sterile pasteur pipettes
Incubator (35 - 37°C), not fan assisted
Refrigerator (2 - 8°C)
Marking Pen
Oxidase strips (MID61G)
Hydrogen Peroxide, use at 3% (w/w), for catalase test see Reference 1
Gram stain reagents
Microscope and Microscope slides
25°C Incubator, not fan assisted

STORAGE

The microwell strips are stable in the unopened foil pouches at 2 - 8°C until the expiry date stated. The Listeria suspending broth and haemolysin reagent should be stored at 2 - 8°C. The haemolysin reagent should be returned to 2 - 8°C immediately after use.

INSTRUCTIONS FOR USE

(Before using this product, refer to Precautions and Limitations)

1. Selection of colonies for identification

- 1.1. Isolates can be tested from any selective or non selective media.
- 1.2. Prior to inoculation into the Microgen Listeria ID, isolates should be checked to ensure they are members of the genus **Listeria**. (short Gram positive bacillus, oxidase negative, catalase positive, motile at 25°C but non motile at 37°C (we recommend that motility be determined by the microscopy method described in Reference 1) Alternatively the Microgen Listeria Latex test (F48) may be employed.

2. Inoculum preparation

- 2.1. Bring the suspending broth to room temperature before inoculation of microwell strips.
- 2.2. Select a single well-isolated colony from an 18-24 hour culture and emulsify it in a vial of Listeria Suspending medium (2.5ml).
- 2.3. Mix thoroughly to produce a homogenous suspension.

3. Inoculation and Incubation

- 3.1 Remove a microwell strip from the foil pouch, place it in the holding frame and remove the lid.
- 3.2 Using a sterile Pasteur pipette transfer 4 drops (approximately 100µl) of the bacterial

suspension to each well of the microwell strip.


- 3.3 As a purity check, transfer 1 drop of the organism suspension onto an appropriate non selective agar plate. Incubate the plate aerobically at 35 - 37°C for 18 - 24 hours.
- 3.4 Add 1 drop of the haemolysin reagent to well 12.
- 3.5 Replace the lid onto the microwell strip and incubate at 35 - 37°C for 18- 24 hours.

4. Interpretation

- 4.1 After incubation remove the lid from the microwell strip and record results on the report forms provided.
- 4.2 Refer to the table of tests (page 1) for guidelines in the interpretation of the results.
- 4.3 The haemolysin reaction should be examined as follows:
 - 4.3.1 Examine the bottom of the well for the presence of a distinct layer or carpet of red blood cells, this should be interpreted as a **NEGATIVE** test. The absence of a distinct layer or carpet of red blood cells should be interpreted as a **POSITIVE** test.
 - 4.3.2 Examine the inoculum in the well. The presence of a clear solution above the distinct layer or carpet of red blood cells in the bottom of the microwell, should be interpreted as a **NEGATIVE** test. The presence of a straw - brown homogeneous solution with the absence of a carpet of red blood cells should be interpreted as a **POSITIVE** test.
 - 4.3.3 If tests are difficult to read, it is suggested that the microwell strip be placed over a light source.



**POSITIVE
HAEMOLYSIS TEST -**
a straw - brown solution
with the absence of a
carpet of
red blood cells



**NEGATIVE
HAEMOLYSIS TEST -**
a clear solution above the
distinct layer or carpet of
red blood cells in the
bottom of the microwell,

- 4.4 Examine the purity plate for viability of the test organism and purity.
- 4.5 The tests on the report form have been organised into triplets (sets of 3 reactions), with each test assigned a numerical value (1,2 or 4). The sum of the positive reactions for each triplet forms a single digit of the Profile Number (Octal Code) that is used to determine the identity of the *Listeria spp.* being identified. The Profile Number (Octal Code) is entered into the Microgen Identification System Software, which generates a report of the five most likely organisms based on the selected database (7).

Report Form

MICROGEN LISTERIA – ID REPORT FORM															
Lab. No. 2894			Specimen Type: GREEN SALAD												
			Date: 28th JANUARY 2002												
	Oxidase	Catalase	Latex Agglut.	Esculin	Mannitol	Xylose	Arabitol	Ribose	Rhamnose	Trehalose	Tagatose	Gluc-1-Phos	M-D-Gluc	M-D-Man	Haemolysis
Reaction															
Result	+	+		+	-	-	+	-	+	+	-	-	+	+	+
Reaction Index	4	2		4	2	1	4	2	1	4	2	1	4	2	1
Sum of Positive Reactions				14			5			4			7		
Profile No: 4547			Final Identification: L. monocytogenes												

WF6031/01/12

PRECAUTIONS

- The Microgen Listeria-ID system is intended for use by qualified laboratory personnel using aseptic technique and appropriate microbiological precautions.
- Used materials must be disposed of safely by autoclaving, incineration or immersion into an appropriate disinfectant prior to disposal.
- The microwell strip lids do not seal the microwells completely so the strips **must not** be incubated in either a CO₂ incubator (due to erroneous pH effects) or fan assisted incubator (potential for excess evaporation)
- The haemolysin reagent should be handled using good microbiological technique to avoid contamination:
 - Always store at 2-8°C
 - Avoid contact of the dropper with the test strip or other surfaces during use and always immediately replace the dropper cap

The haemolysin reagent may not perform properly if it has deteriorated eg. due to heavy contamination in use. Deteriorated reagent should not be used. Signs of deterioration are significant haemolysis of the vial contents or the reagent may appear a dark wine-brown colour.

If the result of the haemolysin test is unclear the isolate should be inoculated on to a sheep blood agar plate and the plate checked for haemolysis after incubation at 35 – 37°C for 18 – 25 hours.

LIMITATIONS

- Although selective media for the isolation of *Listeria spp.* are formulated to inhibit the growth of a wide range of contaminating normal flora, organisms which resemble *Listeria spp.* on these media may grow through (*Bacillus spp.*, *Enterococcus spp.* and *Staphylococcus spp.*).
- The Microgen Listeria ID system has been designed to identify organisms belonging to the genus *Listeria* and no other genera. If the isolate being identified does not hydrolyse Aesculin or ferment Trehalose or Arabitol the gram stain, motility, oxidase and catalase should be re checked.
- Specimens or samples may contain a mixture of species therefore the selection of a single well-isolated colony is critical to obtaining the most accurate result.
- Inoculation of a purity plate is recommended as it will confirm that a single species was inoculated into the test strips.

QUALITY CONTROL

The performance of the Microgen Listeria ID system should be monitored using appropriate control strains. The following are recommended for independent laboratory assessment:

	E S C	M A N	X Y L	A R L	R I B	R H A	T R E	T A G	G I P	M D G	M D M	H E M
<i>L.monocytogenes</i> (ATCC 35152, NCTC 7973)	+	-	-	+	-	+	+	-	-	+	+	+
<i>L.inocua</i> (ATCC 33090, NCTC 11288)	+	-	-	+	-	+	+	-	-	+	+	-
<i>L..grayi</i> (ATCC 19120, NCTC 10815)	+	+	-	+	+	-	+	-	-	-	+	-

DATABASE

	ESC	MAN	XYL	ARL	RIB	RHA	TRE	TAG	GIP	MDG	MDM	HEM
<i>L.monocytogenes</i>	100	0	0	97	0	98	97	0	2	99	98	99
<i>L.inocua</i>	100	0	1	100	0	70	100	0	0	100	100	0
<i>L.welshimeri</i>	100	0	95	100	0	87	100	94	0	98	94	0
<i>L.seeligeri</i>	100	0	100	100	0	0	97	0	0	100	5	93
<i>L.ivanovii</i>	100	0	97	100	42	5	86	0	92	95	0	90
<i>L.grayi</i>	100	97	0	100	100	0	98	0	0	30	94	0

Figures denote percentage positive strains
Highlighted reactions are confirmatory for *Listeria spp.*

REFERENCES

1. On Line Bacteriological Analytical Manual - www.FDA/CFSAN Bacteriological Analytical Manual Online, Chapter 10 - Detection and Enumeration of *Listeria monocytogenes* in Foods.
2. Confirmation of *Listeria* species Method 11.3:1995 CCRFA Microbiological Methods Manual
3. AS/NZS 1766.2.15:1998 Examination for specific organisms – *Listeria monocytogenes* in dairy products.
4. Rodriguez L.D., J.A. Vazquez Boland, j.f. Fernandez Garayzabal, P. Echalecu Tranchant, E. Gomez-Lucia, E.F. Rodriguez Ferri and G. Suarez Fernandez. 1986 A Microplate Technique to Determine Hemolytic Activity for Routine Typing of *Listeria* Strains. **24**:99 – 103.
5. Mira-Gutierrez J. and C.Perz De Lara and M.A. Rodriguez-Igesias. 1990. Identification of species of the genus *Listeria* by fermentation of carbohydrates and enzymatic patterns. Acta Microbiologica Hungarica **37**:123 – 129.
6. Wilkinson B.J. and D.Jones. 1977. A Numerical Taxonomic Survey of *Listeria* and Related Bacteria. J.Gen. Microbiol. **98**: 399 – 421.
7. Lapage S.P, S.Bascombe, W.R. Willcox and M.A.Curtis. 1973 Identification of Bacteria by Computer: General Apects and Perspectives J.Gen. Microbiol. **77**: 273 -290

Microgen Bioproducts Ltd
1 Admiralty Way
Camberley
Surrey
GU15 3DT
U.K.

WF5048/09/2004
