## CHROMagar™ ESBL

Click below:

EN

FR

ES

# ● CHROMagar™ ESBL

**For overnight detection of Gram-negative bacteria  
producing Extended Spectrum Beta-Lactamase**

## For overnight detection of Gram-negative bacteria producing Extended Spectrum Beta-Lactamase

### Background

ESBL (Extended Spectrum  $\beta$ -Lactamases) are enzymes that mediate resistance to penicillins, extended-spectrum third generation cephalosporins (C3G) and monobactams. ESBL-producing *Enterobacteriaceae* started to appear in the 1980s, and have since emerged as some of the most significant hospital-acquired infections with *Escherichia coli* and *Klebsiella spp.* being the main actors, but other Gram-negative species have also been observed. Emergence of ESBL-producing isolates has important clinical and therapeutic implications:

- Resistance determinants for ESBL production are carried on plasmids that can be easily spread from organism to organism.
- The spread of resistance toward extended-spectrum cephalosporins may lead to increased prescription of more broad-spectrum and expensive drugs.
- These resistant isolates may escape detection with routine susceptibility testing performed by a clinical microbiology laboratory, which can result in adverse therapeutic outcomes.

Therefore, the early detection of ESBL-producing bacteria carriers is important to minimise their impact and the spread of infections and customise therapeutic patient treatment.

### Intended Use

CHROMagar™ ESBL is a selective and differential chromogenic culture medium, intended for use in the qualitative direct detection of gastrointestinal colonization with Extended-spectrum-beta-lactamase-resistant Enterobacteria (ESBL) to aid in the prevention and control of ESBL in healthcare settings. The test is performed with rectal swab and stools from patients to screen for ESBL colonization. Results can be interpreted after 18-24 h of aerobic incubation at 35-37 °C. CHROMagar™ ESBL is not intended to diagnose infection caused by ESBL producers nor to guide nor monitor treatment for infections. A lack of growth or the absence of colonies on CHROMagar™ ESBL does not preclude the presence of ESBL-producing bacteria. Further identification, susceptibility testing, and epidemiological typing is needed on suspect colonies.

### Medium Performance


- 1 FAST RESULTS**  
Detection after overnight incubation
- 2 SPECIES DIFFERENTIATION**  
thanks to the chromogenic performances of supplemented CHROMagar™ Orientation. Indeed, the product is composed of a powder base CHROMagar™ Orientation and a supplement to enhance ESBL-producing bacteria.
- 3 HIGH SENSITIVITY AND SPECIFICITY**

	Analytical data *		Clinical data **
	CHROMagar™ ESBL	CHROMagar™ ESBL	Competitor's medium
Sensitivity	100 %	98.3 %	97.5 %
Specificity	97 %	72.3 %	72.9 %

\* Data obtained after a 18-24 h incubation at 37 °C in aerobic conditions in the study « Rapid detection of ESBL-producing Enterobacteria on chromogenic medium: Colorex Orientation/ESBL ». Laudat et al. Poster SFM 2010.

\*\* Data obtained by testing 2 337 rectal screening samples, on plates incubated at 37 °C for 18 h, in the study « Comparison of five media for detection of Extended-Spectrum beta-Lactamase by use of the wrap instrument for automated specimen processing ». Grohs et al., 2013. *J.Clin. Microbiol.*

- 4 TIME AND WORKLOAD SAVINGS**  
Direct culture from specimen. There is no need of a selective pre-enrichment.
- 5 FLEXIBILITY**  
CHROMagar™ ESBL is supplied with a more than 18 months shelf-life . This allows flexibility of use, whether in an epidemic situation with many patients to screen, or whether for random surveillance of cultures.

### Plate Reading

- *E. coli* ESBL  
→ dark pink to reddish
- *Klebsiella, Enterobacter, Citrobacter* ESBL  
→ metallic blue (+/- red halo)
- *Proteus* ESBL  
→ brown halo
- *Pseudomonas* ESBL  
→ translucent cream to blue
- *Acinetobacter* ESBL  
→ cream, opaque

### Medium Description

<b>Powder Base</b> (CHROMagar™ Orientation)	Total ..... 33 g/L Agar ..... 15.0 Peptone and yeast extract ..... 17.0 Chromogenic mix ..... 1.0 Storage at 15/30 °C - pH: 7.0 +/-0.2 Shelf Life ..... > 18 months
<b>+</b>	
<b>CHROMagar™ ESBL supplement</b> (included in the pack)	Selective mix (Powder form) ..... 0.57 g/L Storage at 2/8 °C Shelf Life ..... > 18 months

Usual Samples	Rectal swabs and stools
Procedure	Direct Streaking. Incubation 18-24h at 35-37 °C. Aerobic conditions

Scientific Publications on this product: available on [www.CHROMagar.com](http://www.CHROMagar.com)  
Please read carefully the instructions for use (IFU document) available on [www.CHROMagar.com](http://www.CHROMagar.com)

### Order References

Please use these product references when contacting your local distributor:


- 5000 mL pack ..... ESRT2 (included in this reference: powder base RT412 + supplement ES372)
- 25 L pack ..... ESRT3-25 (included in this reference: powder base RT413-25 + supplement ES373-25)

Manufacturer: CHROMagar  
4 place du 18 juin 1940 75006 Paris - France  
Email: [CHROMagar@CHROMagar.com](mailto:CHROMagar@CHROMagar.com)  
Website: [www.CHROMagar.com](http://www.CHROMagar.com)  
Find your nearest distributor on  
[www.CHROMagar.com/contact](http://www.CHROMagar.com/contact)

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