


CHROMagar™ Vibrio

Save time, save labour costs
with the original
CHROMagar



CHROMagar Paris France - Fax : +331 45 48 06 06
CHROMagar™ is a trademark of Dr. A. Rambach
LF-EXT-008F-Feb 2005

The Chromogenic Media Pioneer

CHR  Magar

www.chromagar.com



CHROMagar™ Vibrio

A *Vibrio* medium with colour differentiation
V.parahaemolyticus, *V.vulnificus*, *V.cholerae*

SPECIES	COLONY COLOUR
<i>V.parahaemolyticus</i>	Mauve
<i>V.vulnificus</i> & <i>V.cholerae</i>	Blue
<i>V.alginolyticus</i>	Colourless

➔ **CHROMagar™ Vibrio** can permit an improved detection of *V.parahaemolyticus* (Kudo *et al* 2001).

➔ **CHROMagar™ Vibrio** has a colony differentiation by colour (mauve, blue, colourless) **more powerful than classical TCBS medium.**

Reference: Yukiko Hara-Kudo *et al.* 2001. Improved Method for Detection of *Vibrio parahaemolyticus* in Seafood. Applied & Environmental Microbiology. December 2001, Vol.67, N°12 p. 5819-5823.

CHROMagar™ Candida	for differentiation of <i>C.albicans</i> , <i>C.tropicalis</i> , <i>C.krusei</i>
CHROMagar™ Orientation	for differentiation of urinary tract pathogens
CHROMagar™ E.coli	for detection and enumeration of <i>E.coli</i>
CHROMagar™ ECC	for detection and enumeration of <i>E.coli</i> and coliforms
CHROMagar™ Liquid ECC	broth for pad technique for <i>E.coli</i> -coliforms (water samples)
CHROMagar™ O157	for detection of <i>E.coli</i> O157
CHROMagar™ Salmonella	for detection of <i>Salmonella</i> including <i>S.typhi</i>
Rambach™ Agar	for detection of <i>Salmonella</i> spp
CHROMagar™ Staph aureus	for detection of <i>Staphylococcus aureus</i>
CHROMagar™ MRSA	for detection of Methicillin Resistant <i>Staphylococcus aureus</i> including low level MRSA
CHROMagar™ Listeria	for detection and enumeration of <i>Listeria monocytogenes</i>
CHROMagar™ Vibrio	for detection of <i>Vibrio parahaemolyticus</i> , <i>V.vulnificus</i> , <i>V.cholerae</i>
CHROMagar™ Pseudomonas	for detection of <i>Pseudomonas</i>

Visit CHROMagar on internet via chromagar.com CHROMagar™ and Rambach™ are trademarks of Dr. A. Rambach

CHROMagar
 4, place du 18 Juin 1940
 75006 Paris – France
 Tel: +33 1 45 48 07 07
 Fax: +33 1 45 48 06 06

The Pioneer of Chromogenic Media

CHR  **Magar™**
 Microbiology