

Rock-IT & Roll

Laboratory Shaker

Complete extraction in a fraction of the time

工欲善其事必先利其器 *Shaker*

可調式自動關機時間
桌上式使用方便
一次可容納20組檢體

除去手動式的疲勞
萃取處理最適用







VOLUME PRICING AVAILABLE UPON REQUEST



800/234-5333 or 517/372-9200
620 Leshar Pl., Lansing, MI 48912
foodsafety@neogen.com • www.neogen.com





Rockit Evaluation Conclusions

Aflatoxin

1 minute Rockit shake is equivalent to a 2 minute blend extraction

Fumonisin

1 minute shake is equivalent to a 2 minute blend extraction

Deoxynivalenol

30 second Rockit shake is equivalent to a 3 minute hand shake

Benefits

- RAPID EXTRACTION TIME – 30 SECONDS TO 1 MINUTE
- EXTRACT 12 SAMPLES AT ONCE
- UTILIZES DISPOSABLE EXTRACTION CONTAINERS
- INEXPENSIVE COMPACT DESIGN FOR BENCH TOP USE
- ELIMINATES HAND SHAKING USED FOR ELISA KIT EXTRACTIONS
- SYRINGE FILTER RAPIDLY CLARIFIES SAMPLE EXTRACTS PRIOR TO ANALYSIS

經美國密蘇里州Trilogy分析實驗室實驗證實

NEOGEN Rock it Roll shaker振盪式萃取比手動搖晃式速度快6倍 (DON)
 NEOGEN Rock it Roll shaker振盪式萃取比均質機速度快2倍 (Aflatoxin)
 NEOGEN Rock it Roll shaker振盪式萃取比均質機速度快2倍 (Fumonisin)