

**We Stand Behind Our Results**

*Choose Neogen.*

## **NeoColumn for Aflatoxin DR Test System**

**精確黃麴毒素免疫親和管**

### **Introducing Neogen's NEW NeoColumn for Aflatoxin DR and DRF 2100 Fluorometer**

NeoColumn for Aflatoxin DR is a direct read column with the same easy steps used with other direct read methods. NeoColumn run samples can be interpreted with Neogen's DRF 2100 fluorometer (or other similar calibrated fluorometers) or can be used as a clean-up for HPLC analysis. The NeoColumn direct read method is USDA-GIPSA approved and has proven accuracy compared to HPLC and other fluorometric methods.



**高效率免疫親和管**

**NeoColumn for Aflatoxin DR**

Direct read (DR), USDA GIPSA approved

Easy to use single use immunoaffinity column for the detection of aflatoxin B<sub>1</sub>, B<sub>2</sub>, G<sub>1</sub>, and G<sub>2</sub>.

### **Direct Read Fluorometer (DRF 2100)**

The easiest, most intuitive fluorometer available.

**高效率螢光光度計**



To order, call 800/234-5333 or visit [www.neogen.com](http://www.neogen.com)

**ISO**  
certified  
9001:2000

**NEOGEN總代理 TPM台灣現用股份有限公司**

[www.tpm-tpm.com.tw](http://www.tpm-tpm.com.tw) TEL: 02-2796 8854 FAX: 02-2796 8039