


# 光學讀取總生菌及大腸桿菌計數儀



Holds 128 sample vials in four incubators.  
Capacity of 128 tests at 4 independent temps (each incubator holds 32 vials). Easy to insert and remove vials.



## Soleris 32

## Soleris 128

### 11種測試瓶功能最強

1. 總生菌
2. 大腸桿菌群
3. 大腸桿菌
4. 腸內菌
5. 酵母菌
6. 乳酸菌
7. 真菌
8. 綠膿桿菌
9. 李斯特菌
10. 金黃色葡萄球菌
11. 無菌測試瓶



總生菌檢測



真菌檢測

### 1024 檢測數量最大

#### Soleris 32

Height: 16 cm  
Width: 45 cm  
Depth: 55 cm  
Weight: 12 kg


#### Soleris 128

Height: 42 cm  
Width: 73 cm  
Depth: 69 cm  
Weight: 60 kg


### 8 台機器連線器



### 2007最新電腦軟體




### 菌落CFU 自動計算



Enter Spec. Level (CFU): 1000

Cutoff Time: 11

### 檢測時間歷史統計圖



Product: beef

Test: beef cc

From: May 2004

### 歷史趨勢分析圖


Start Individual Samples

Prepare Next Run | Start Prepared Run

Start Individual Samples | Stop Individual Samples | Archive Individual Curves

1A1	1	2	3	4
A	Stopped	Stopped	Stopped	Stopped
B	Stopped	Stopped	Stopped	Stopped
C	Stopped	Stopped	Stopped	Stopped
D	Stopped	Stopped	Stopped	Stopped
E	Stopped	Stopped	Stopped	Stopped
F	Stopped	Stopped	Stopped	Stopped
G	Stopped	Stopped	Stopped	Stopped
H	Stopped	Stopped	Stopped	Stopped

Start From : [ ] To : [ ] [Undo] [OK]



# SOLERIS™

## 強大軟體 管理系統



**超強偵測力 每5分鐘檢測一次  
不用數菌落 6-12小時直接發報告**



李斯特菌檢測瓶

Graph From Archive [Exit]

1A1	1	2	3	4
A	26Apr1A1 Ground beef GB-TVC	03AprMD2 Milk Milk TVC 6.9 2.1E+4	03AprMD3 Milk Milk TVC 9.0 1.4E+3	03AprMD4 Milk Milk TVC 7.7 7.9E+3
B	03AprMD1 Milk Milk Coliform 5.9 1.7E+5	03AprMD2 Milk Milk Coliform 10.7 1.6E+2	03AprMD3 Milk Milk Coliform 9.6 7.8E+2	03AprMD4 Milk Milk Coliform
C	03AprYD1 Yogurt Yogurt Coliform	03AprYD2 Yogurt Yogurt Coliform 5.0	03AprYD3 Yogurt Yogurt Coliform 9.4	03AprYD4 Yogurt Yogurt Coliform
D	03AprYD1 Yogurt yeast	03AprYD2 Yogurt yeast	03AprYD3 Yogurt yeast 12.0	03AprYD4 Yogurt yeast
E	03AprPD1 Pork Pork-TVC 5.1 5.0E+4	03AprPD2 Pork Pork-TVC 7.5 1.1E+6	03AprPD3 Pork Pork-TVC 9.1 3.7E+5	03AprPD4 Pork Pork-TVC 11.7 6.1E+4
F	03AprPD1 Pork Pork Coliform 6.8 2.7E+5	03AprPD2 Pork Pork Coliform 8.9 3.4E+4	03AprPD3 Pork Pork Coliform 11.0 4.2E+3	03AprPD4 Pork Pork Coliform
G	03AprCD1 ChickenRinse ChickenRinse TVC 9.2 5.9E+4	03AprCD2 ChickenRinse ChickenRinse TVC 4.5 7.2E+6	03AprCD3 ChickenRinse ChickenRinse TVC 5.6 2.3E+6	03AprCD4 ChickenRinse ChickenRinse TVC 7.2 4.6E+5
H	03AprCD1 ChickenRinse ChickenRinse Colifo 5.1 4.6E+5	03AprCD2 ChickenRinse ChickenRinse Colifo 7.7 6.3E+4	03AprCD3 ChickenRinse ChickenRinse Colifo 9.7 6.4E+3	26Apr1H4 Past. Milk PMilk-TVC

Chosen Locations: [ ] [Undo] [Plot]

