

USP CHAPTERS 61 & 62 *effective May 2009*

USP/EP/JP Compliance Performance testing for all Acumedia products have been revised to meet or exceed all USP harmonized specifications

Test Medium	Property	Test Strains	Cat.	Type
Test for Bile-Tolerant Gram Negative Bacteria				
Enterobacteriaceae Enrichment Broth, Mossel	Growth	<i>E. coli</i>	7603	Dehydrated 500g
		<i>P. aeruginosa</i>		
	Inhibitory	<i>S. aureus</i>		
Violet Red Bile Glucose (VRBG) Agar	Growth	<i>E. coli</i>	7425	Dehydrated 500g
		<i>P. aeruginosa</i>		
Test for Escherichia coli				
MacConkey Broth	Growth	<i>E. coli</i>	7185	Dehydrated 500g
	Inhibitory	<i>S. aureus</i>		
MacConkey Agar	Growth	<i>E. coli</i>	7102	Dehydrated 500g
	Inhibitory	<i>S. aureus</i>		
Test for Salmonella				
Raooaort Vassiliadis Salmonella Enrichment Broth	Growth	<i>S. enterica serotype</i>	7730	Dehydrated 500g
	Inhibitory	<i>S. aureus</i>		
Xylose-Lysine-Deoxycholate Agar (XLD)	Growth	<i>S. enterica serotype</i>	7166	Dehydrated 500g
	Indicative	<i>E. coli</i>		
Test for Pseudomonas aeruginosa				
Cetrimide Agar	Growth	<i>P. aeruginosa</i>	7688	Dehydrated 500g
	Inhibitory	<i>E. coli</i>		
Test for Staphylococcus aureus				
Mannitol Salt Agar (MSA)	Growth	<i>S. aureus</i>	7143	Dehydrated 500g
	Inhibitory	<i>E. coli</i>		
Test for Clostridia				
Reinforced Medium for Clostridia	Growth	<i>Clostridium sporogenes</i>	7714	Dehydrated 500g
Columbia Agar	Growth	<i>Clostridium sporogenes</i>	7125	Dehydrated 500g
Test for Candida albicans				
Sabouraud Dextrose Broth	Growth	<i>C. albicans</i>	7617	Dehydrated 500g
		<i>A. brasiliensis</i>		
Sabouraud Dextrose Agar	Growth	<i>C. albicans</i>	7150	Dehydrated 500g
		<i>A. brasiliensis</i>		
Microbial Enumeration Tests		Type	<h1>USP</h1> <h2>Harmonization</h2>	
Tryptic Soy Broth (TSB)		Dehydrated 500g		
Tryptic Soy Agar (TSA)		Dehydrated 500g		
Potato Dextrose Agar (PDA)		Dehydrated 500g		

Acumedia Meeting Regulatory Requirements

TPM 台灣現用股份有限公司 www.tpm-tpm.com.tw TEL: 2796 8854 2796 2025 Fax: 2796 8039

Acumedia 台灣總代理 e-mail: tpm.tpm@msa.hinet.net